Dentimplant

Abutment Surgical Use Guide Abutment

- 1- Take an impression after soft tissue is healed.
- 2- Choose the right type and size of abutment by comparing the radiography with the taken impression and by considering adjusting teeth.
- 3- Connect the chosen abutment to the fixture in the mouth.
- 4- Take an impression at either fixture level or at abutment level in order to produce the fixture prosthesis.
- 5- In the case of cement-retained prosthesis, use a hard finishing plaster to do a wax-up on the abutment, and manufacture the prosthesis.
- 6- Attach the final prosthesis to the implant or the abutment in the mouth and connect it to the fixture and then check the occlusion between the upper and lower jaws.

Usage Guidelines

1. Patient Population

Dental implants are designed for use in patients with complete or partial edentulism, who have completed their growth and development, do not have the conditions specified in the contraindications. (Patients older than 18 years of age who have completed awbone development.)

2. Indications

Dentcon abutment and prosthetic parts are used in lower and/or upper jaw to support crown, bridge or total prostheses in treatment of full or partial missing of teeth.

3. Contraindications

- 1- Patients with a history of heart failure or artery hardening;
- 2- Uncooperative patients or patients with mental or physical disorder which may result in unstable implant, failure to fix the implant or complications in management after operation;
- 3- Patients with the problem in bone remodeling, microcirculation or blood clotting function;
- 4- Pregnant patients;
- 5- Patients with high blood pressure or diabetes mellitus;
- 6- Patients with bone disease (osteoporosis, osteomalacia) or bone metabolism disorders;
- 7- Patients allergic to Titanium

4. Pay Attention

- 1.It should be done by specialist physicians because bone tissue or transplant may be damaged due to the wrong surgical method.
- 2.This product should not be reused and should be applied for its original purpose.
- 3.Misused or damaged product should be discarded.
- 4.If the product is faulty, it must be sent back.
- 5.It should be used carefully to avoid damage and deformation of the product.

5. Side Effect

Swelling , Local inflammations , Gingival injuries, Local pain, Micro hemorrhages, Aspiration or swallowing of parts during operation, Triggering of the pharyngeal reflex during the insertion of the product, Screw fracture, Screw loosening, Component or prosthesis failure/fracture, Hypersensitivity and/or allergic reactions, Other toxicity reactions

6. Warnings

This product is not sterile and must be sterilized before placing in the patient's mouth. When the dentist places the final denture in the patient's mouth, he or she should check the osseointegration with x-rays, a stability checker, and percussion tests. Abutment screws should be tightened to the specified torque values. If the screw is tightened at a value within this torque range, this may damage the screw and its seat, causing it to break. Tightening at a low torque can also damage it. When connecting the screw, be especially careful not to use too much force until you are sure that the screw is seated correctly. If the initial placement of the screw is improper, it may result in damage to the screw, the abutment, or the entire implant structure. The chewing force of the prosthesis should not be applied to one side. Otherwise, stress fractures may occur.

7. Frequency of Occurrence of the Most Common Identified Risks

- Implant Loss / Osseointegration Failure:

Implant survival is typically 96%–99%.
Implant loss may occur at 0%–3% in routine cases and 8%–22% in complex cases.

- Peri-implantitis / Infection:

Peri-implant infection rates are generally 0%–2.7% and are associated with oral hygiene.

- Marginal Bone Loss:

Marginal bone loss is typically 0.02–1.22 mm with an average of 0.4–0.7 mm.

- Soft Tissue Recession / Esthetic Complications:

≥1 mm recession may occur in up to 11% of cases; thin biotype is a key risk factor.

- Mechanical / Prosthetic Complications:

Mechanical/prosthetic complications occur at 1%–12% and are manageable with maintenance.

8. Sterilization

Abutments are supplied to the market as non-sterile. Therefore, before abutments are used, sterilize in an autoclave at 121 °C for 15 minutes and wait until it dries.

9. Magnetic Resonance Imaging (MRI) Safety Information

- * Static magnetic fields of 1.5 and 3.0 Tesla only
- * Magnetic spatial dropout field of up to 970 Gauss/cm (9.7 T/m) or less
- * Whole body average specific absorption rate (SAR) of 2 W/kg (Normal Use Mode)
- * Under the scanning conditions specified above, our products are expected to produce a maximum temperature increase of 3°C after 15 minutes of continuous scanning.

10. More Information

In case of product complaints or for more information please contact.

Summary of Safety and Clinical Performance (SSCP) will be available in the European database. This database can be accessed at:
<https://ec.europa.eu/tools/eudamed> The SSCP document related to the DENTIMPLANT abutment/implant system can be accessed through this database. Reference document: SSCP01 Dental Implant System Summary of Safety and Clinical Performance

NOTE: You can view the electronic version of the manuals at [www.dentimplant.com.tr] (<http://www.dentimplant.com.tr>)

To access the SSCP file:
<https://www.dentimplant.com.tr/files/>

11. Manufacturer Information


Company Name: DENTIMPLANT
MEDİKAL SAĞLIK İÇ DİŞ TİCARET SAN. A.Ş.

Company Address: Ostim OSB Mah. 1269. Cad. No:17
Yenimahalle / ANKARA

12. Contact

If any adverse effects, complaints occur, pls contact us with the following informations as below:
Contact Person PRRC : Seçil KOCABIYIK
0312 424 10 28

Website: www.dentimplant.com.tr

 Manufacturer Identification	 Production Date
 Expiry Date	 Lot number
 Reference Number	 Non Sterile
 Do not Resterilize	 Do not Use if Package Is Damaged
 Keep Away from the Sun	 For Single Use Only
 Caution	 Check the User Manual
 MR Conditional	 CE 2696 Conformity Mark for Medical Device Regulation 2017/745
 Steam Sterilization	

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